



Applications

- Weighing
- Parts counting
- Percentage weighing
- Check weighing
- Density determination
- Below balance weighing

Features

- 9 weighing units
- Date and time
- Full range tare
- Zero Tracking
- Multilingual display
- Large backlit display with dual text prompts
- Capacity tracker
- Large stainless steel pan
- Solid metal housing
- Overload protection
- Lock down mounting slot for Kensington™ type lock and cable
- Non-slip adjustable levelling feet
- Colour coded, sealed keypad
- Dual tare keys
- Removable draft shield supplied on 0.001g readabilities
- Modern low profile design
- RS-232 bi-directional interface
- GLP print outs
- External calibration
- Rechargeable battery
- AC Adapter
- Auto sleep / power down function to save battery life

Technical Specifications

Model	PGL 20001
Maximum capacity	20kg
Readability	0.1g
Repeatability (S.D.)	0.2g
Linearity (±)	0.4g
Pan size	400×300mm

General Specification

Typical Stabilisation time	2-3 seconds
Units of Measurement	Gram (g), Kilogram (kg), Carat (ct), Grain (GN), Newton (N), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz), Custom unit
Operating Temp	10°C 40°C
Power supply	AC adaptor and rechargeable battery pack
Calibration	External
Display	Backlit LCD with dual digits and capacity tracker
Overall Dimensions	358x251x104mm
Net Weight	6kg

Accessories

- No.8023 Adam Printer
- No.9061 AdamDU - Data Collection Program
- No.9028 RS-232 Cable
- No.8036 Anti-Vibration Table
- No.7995 Below Balance Hanger
- No.8031 In-use Cover
- No.8002 Dust Cover
- No.8030 Security Lock and Cable
- No.3584 Thermal Paper Roll for Adam Printer
- No.7453 Density Determination Kit (mg models only)
- No.F1008 F1 200g Calibration Weight for PGL 203
- No.F1008 F1 200g Calibration Weight for PGL 303
- No.F1010 F1 1kg Calibration Weight for PGL 2002
- No.F1010 F1 1kg Calibration Weight for PGL 3002
- No.F1011 F1 2kg Calibration Weight for PGL 4001
- No.F1012 F1 5kg Calibration Weight for PGL 6001
- No.F1012 F1 5kg Calibration Weight for PGL 8001
- No.F1012 F1 5kg Calibration Weight for PGL 10001
- No.F1012 2 x F1 5kg Calibration Weights for PGL 20001